

Autoradiography Film

RayMax™ Autoradiography Film

RayMax™ Beta

Introducing the latest addition to the ICN family of radiochemicals and products for molecular biology research – RayMax™ autoradiography film!

RayMax™ Universal – a double coated, multipurpose film for general autoradiography and fluorography offering high resolution results at economical prices.

The film consists of a clear plastic base coated on both sides with a photographic emulsion. The emulsion is coated with a protective anti-scratch layer.

RayMax™ Tritium – a specialized single coated film without an anti-scratch layer for direct autoradiography of tritium or ¹²⁵I.

RayMax™ Beta – a single coated, high silver content film for high resolution direct autoradiography of ¹²⁵I, ³³P, ³⁵S or ¹⁴C.

RayMax™ Universal

Cat. No.	Film Size	Qty.	Price
821716	18 x 24 cm	25 sheets	97.20
821717	18 x 43 cm	25 sheets	132.30
821718	24 x 30 cm	25 sheets	139.50
821719	30 x 40 cm	25 sheets	171.90
821720	35 x 43 cm	25 sheets	213.30

RayMax™ Tritium

Cat. No.	Film Size	Qty.	Price
821711	18 x 24 cm	10 sheets	125.10
821715	24 x 30 cm	10 sheets	204.30

Cat. No.	Film Size	Qty.	Price
821722	18 x 24 cm	25 Sheets	52.20
821723	18 x 43 cm	25 Sheets	127.80
821724	24 x 30 cm	25 Sheets	67.50
821725	30 x 40 cm	25 Sheets	146.70
821726	35 x 43 cm	25 Sheets	151.20
821727	5" x 7"	25 Sheets	42.30
821728	8" x 10"	25 Sheets	52.20
821729	18 x 24 cm	75 Sheets	142.20
821730	35 x 43 cm	75 Sheets	382.50
821731	8" x 10"	75 Sheets	142.20
821732	18 x 24 cm	50 Envelopes	78.30
821733	35 x 43 cm	50 Envelopes	292.50
821734	8" x 10"	50 Envelopes	117.00